



The GLOVE™ Foot is a great 'add-on' or stand alone additional leg support during fitting of sling or when applying dressings.

Easily slipped over heel with simple strap attachment direct to hoist spreader bar, the GLOVE™ Foot provides improved comfort to the thigh area and reduces excessive load for users in a leg cast. It can also help reduce back injury risk to carers during routine care giving or whilst aiding into/out of bed as well as useful during therapy treatments.

FEATURES & BENEFITS

- Reduce leg load for improved comfort
- Helps reduce carer back injury risk
- Highly supportive
- Easy to attach and remove
- Adjustable strap for multi-positioning of leg
- Offered in three sizes
- Compatible with all sling types





Specialist SLING STYLE PRIMARY USE Leg support

TRANSFER

Paediatric / Adult / Bariatric SIZE SWL (KG) 75

ATTACHMENT TYPE Adjustable strap

CLIP OPTION N/A MODESTY COVER

PICK UP Underfoot POSITIONING AID N/A

HEAD SUPPORT N/A APERTURE N/A AMPUTEE

Soft-knit Polyester MATERIALS

×

PADDING

STANDARD COLOUR

ACCESSORIES × CUSTOMISABLE

Lifetime WARRANTY

COMPATIBILITY All glove sling™ variants OTHER INFO Supplied as single unit

Call 033 012 80808 to book a free assessment